Chirurgische Instrumente

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General Surgical Instruments

Instructions for Use

Instrumentos quirúrgicos generales

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Surgical Instruments – Instructions for use

1. General Notices

Before using the instruments for the first time, as well as before each subsequent use and before returning them for repair or servicing, they must be cleaned, disinfected and sterilized in accordance with our instrument processing instructions, which are available on the Internet at www.avellanus.com but can also be obtained from Avellanus on request.

2. Intended Purpose

The instruments may only be used for their intended purpose in the specified medical fields, with use being restricted to adequately trained and qualified personnel. The treating physician or user is responsible for selecting the right instrument(s) for the surgical task/application at hand as well as for their safe handling. This includes ensuring an adequate level of training, knowledge and experience.

3. Inspection & Testing Prior to Reuse

Before each use, the instruments must be thoroughly inspected for damage such as fractures, cracks or deformation, as well as for functional reliability. Special attention must be paid to cutting edges, tips, joints, box locks, ratchets and all movable parts. If wear, corrosion, deformation, porosity or other damage is detected, the instrument must be immediately withdrawn from service.

4. Processing

Due to their alloy, stainless steel instruments typically develop a passive film in the form of a protective layer. However, this film does not protect them well against chemical attack by chloride ions and aggressive media and liquids! Therefore, in addition to the instrument manufacturer's endeavors to select the right materials and process them carefully, the user must make an important contribution by ensuring proper instrument processing along with adequate and regular care.

To treat our reusable surgical instruments correctly, we recommend using the following procedures:

4.1. Machine Cleaning

It is always preferable to use machines (washer-disinfectors) for instrument cleaning and disinfection because, unlike manual procedures, machine processes can be easily standardized.

Be sure to observe and follow the operating and loading instructions provided by the machine manufacturer. In addition, only the cleaning agents recommended by the manufacturer should be used for the application at hand.

- Jointed instruments should always be processed in open condition. Be sure to arrange the items so that the water can easily flow out of cannulations, blind holes and cavities.
- Complex instruments must be totally taken apart before cleaning.
- For instruments with long or narrow lumens, standard procedures should be used only if the hot disinfectant can easily flow through the lumens and safe rinsing is guaranteed.
- The instrument trays used for cleaning must always be loaded correctly to ensure proper cleaning.
- Store the instruments correctly in the tray. Be sure to prevent "rinsing shadows"!
- For ergonomic reasons, a total tray weight of 8–10 kg should not be exceeded.
- When removing the instruments from the machine after cleaning, be sure to check them for cleanliness (visible dirt). This especially applies to cannulated instruments or those with blind holes. If necessary, repeat the cleaning cycle or pre-clean manually as specified in section 4.2 – Manual Cleaning.
- The final rinse must guarantee that any residues left after the cleaning stage are completely removed. Be sure to use only fully demineralized or distilled water for the final rinse!

4.2. Manual Cleaning

It is best to disinfect and clean the instruments immediately after use in accordance with our recommendations for surgical instrument processing and sterilization. This avoids instrument incrustations that would make cleaning and disinfecting more difficult.

The following should be observed in manual cleaning processes:

- The solutions used for manual cleaning must always be prepared in accordance with the manufacturer's instructions.
- Use a suitable brush for cleaning lumens, cannulations, blind holes and cavities, making sure that every part of the inner surface can be properly accessed.
- Use a soft brush and a neutral or mildly alkaline detergent for removing blood and other residues.
- Never use metal brushes or metal sponges for manual cleaning.
- To ensure proper instrument functioning, verify that all movable parts have been thoroughly cleaned.
- Clean jointed instruments in closed as well as open condition.
- Take instruments fully apart where applicable.
- Pay special attention to slots, ratchets, joints and box locks, narrow lumens, blind holes and other areas that are hard to access.
- Suitable trays or baskets (e.g. sterilization trays or wire baskets) must be used for storing surgical instruments properly during the cleaning process.

4.3. Ultrasonic Treatment

Effective ultrasonic cleaning requires placing the surgical instruments into suitable sterilization trays or wire baskets in an open condition. As the use of plain warm water alone cannot produce satisfactory results, it is necessary to add a suitable cleaning agent. Be sure to observe and follow the manufacturer's instructions with regard to concentration as well as temperature.

Since an excessive dirt content of the cleaning solution has an adverse effect on the cleaning results, the solution must be replaced at regular intervals in accordance with the manufacturer's instructions. In the same manner, prescribed immersion or ultrasonic treatment times must be strictly observed. As a rule, ultrasonic cleaning must always be followed by a rinsing cycle. Be sure to check the instruments (where applicable) for loosened components after the ultrasonic bath. To prevent water spots ("spotting"), fully demineralized or distilled water should be used for the final rinse.

4.4. Chemical Disinfection

1. The solutions employed for chemical disinfection must always be used in accordance with the manufacturer's instructions.
2. Pure water must be used for preparing the dilutions specified for the chemical disinfectants. The addition of detergents is not permitted. The manufacturer's instructions regarding exposure times and concentration must be duly observed in each case.
3. Disinfecting solutions must be prepared afresh every day. Extended/multiple use can easily lead to the following problems:
 - a. Increased concentration due to evaporation (corrosion risk)
 - b. Excessive dirt load (corrosion risk plus lower effectiveness)
4. Following disinfection, it is important to rinse all items sufficiently under clear running water. To prevent the formation of water spots, the use of fully demineralized water is highly recommended.
5. Immediately after completion of the cleaning and rinsing cycles, surgical instruments must be sufficiently dried.

5. Inspection After Cleaning

Following cleaning, the instruments must be macroscopically clean, i.e. free from visible dirt or deposits. Instruments with stains or spots must be withdrawn from service at once and given special treatment. All movable parts, working tips and (scissor) blades should be inspected with particular care.

If damage or malfunction is detected, the instrument must also be withdrawn from service immediately.

6. Instrument Care

The term "care" refers to the application of instrument oil/milk to the surfaces of surgical instruments (especially to movable parts such as joints). "Instrument oil" means physiologically safe paraffin oil as specified by the German Pharmacopoeia (DAB 8) or equivalent regulations (European or US Pharmacopoeia), while "instrument milk" refers to a white-oil-in-water emulsion.

As a rule, surgical instruments must be subjected to regular care, which means each time before a functional test is carried out. At the same time, it is important to prevent "gumming" of the joints due to an accumulative effect, especially in instruments that are continuously in use.

7. Sterilization Packaging

Prior to sterilization, the instruments must be adequately packaged, e.g. using containers that meet DIN EN 868-8 requirements, or wrapping them with suitable packaging materials.

The packaging method used must comply with the relevant standards (e.g. the EN 868 or ISO 11607 series of standards).

8. Steam Sterilization

Sterilization must be performed according to a validated steam sterilization method, e.g. using a sterilizer satisfying EN 285/ANSI/AAMI/ISO 11134-1993 or ANSI/AAMI ST46-1993 requirements, validated in accordance with EN 554/ISO 13683.

It is essential to keep the steam free from any foreign bodies such as rust particles and other impurities. This helps to prevent instrument corrosion or surface damage caused by deposits.

The steam used for sterilization must comply with DIN 58946, Part 7. The user instructions provided by the steam sterilizer manufacturer must be duly observed. Instruments incorporating locks or ratchets must be sterilized in an open condition or with the ratchet set to the first notch.

In specific, the following sources should be observed for cleaning, disinfection and sterilization:

- DIN EN ISO 17664: "Sterilization of medical devices – Information to be provided by the manufacturer for the processing of resterilizable medical devices"
- EN 285: "Sterilization – Steam sterilizers – Large sterilizers"
- EN 554: "Sterilization of medical devices – Validation and routine control of sterilization by moist heat"
- DIN EN 556-1: "Sterilization of medical devices – Requirements for medical devices to be designated 'STERILE' – Part 1: Requirements for terminally sterilized medical devices"
- DIN 58946-7:2004-09: "Sterilization – Steam sterilizers – Part 7: Requirements on installation and services"
- "Proper Maintenance of Instruments" (Red Brochure)
- AKI (Arbeitskreis Instrumentenaufbereitung/"Instrument Preparation" Working Group): <http://www.a-k-i.org> – "Empfehlungen für die Validierung und Routineüberwachung von Sterilisationsprozessen mit feuchter Hitze für Medizinprodukte" (Recommendations for the validation and routine control of processes using moist heat for the sterilization of medical devices)
- DGKH recommendation: <http://www.dgkh.de> – "Anforderungen an die Hygiene bei der Aufbereitung von Medizinprodukten" (Hygiene requirements for the processing of medical devices)
- Recommendation of the Hospital Hygiene and Infection Prevention commission of the Robert Koch Institute (RKI) and of the Federal Institute for Drugs and Medical Devices

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9. Processing Restrictions, Disposal

Frequent reprocessing has only a minor influence on the service life of surgical instruments, which is mainly a function of use-related wear and tear or damage. Please dispose of your instruments in accordance with relevant local regulations, or have them properly recycled, once they have reached the end of their life cycle.

10. Disclaimer

Safety notice: The responsibility for proper instrument cleaning, disinfection and sterilization lies with the operator/user of the product. Be sure to observe and follow all relevant national regulations (incl. relevant restrictions)!

Avellanus, as the manufacturer of the products, accepts no liability for direct or consequential damage caused by improper use, handling, processing, sterilization or maintenance. Unauthorized instrument repair (by firms or persons not specifically authorized by Avellanus to perform such work) shall void the warranty given!

Non-observance of these notices, as well as improper handling or use of products supplied by us, will void your rights under the warranty. Consequently, Avellanus shall not be liable for any resulting damage in such cases.

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Avellanus med GmbH