

## ENDOSCOPY



# ARTHROSCOPY PUNCHES

POWERFUL PINLESS HINGE DESIGN

► Top Sellers



STELLA **LIGHT**  
Coated light weighed aluminium/steel handle

STELLA **STANDARD**  
Stainless steel handle



STELLA **LOCK**  
Stainless steel handle with slide lock

# ARTHROSCOPY PUNCHES



	Description	Tip profile in mm (A)	Bite width in mm (B)	Tip width in mm (C)	STELLA STANDARD	STELLA LIGHT
<b>WIDE PUNCH LOW PROFILE</b>						
	Straight medium	1.90	3.20	5.00	578-020-00	577-020-00
	Curved left 12°; medium	1.90	3.20	5.00	578-035-03	577-035-03
	Curved right 12°; medium	1.90	3.20	5.00	578-040-04	577-040-04
	Curved up 15°; medium	1.90	3.20	5.00	578-030-02	577-030-02
<b>WIDE PUNCH LOW PROFILE; UPBITER</b>						
	Upbiter straight; medium	1.90	3.20	5.00	578-025-10	577-025-10
<b>WIDE PUNCH STANDARD PROFILE</b>						
	Straight medium	2.30	3.20	5.00	578-120-00	577-120-00
	Curved up 15°; medium	2.30	3.20	5.00	578-130-02	577-130-02
<b>WIDE PUNCH STANDARD PROFILE; UPBITER</b>						
	Upbiter straight; medium	2.30	3.20	5.00	578-125-10	577-125-10
	Upbiter curved left 12°; medium	2.30	3.20	5.00	578-135-13	577-135-13
	Upbiter curved right 12°; medium	2.30	3.20	5.00	578-140-14	577-140-14
<b>BASKET PUNCH SLIM</b>						
	Straight	1.90	1.70	3.10	578-210-00	577-210-00
	Curved up 15°	1.90	1.70	3.10	578-215-10	577-215-10
<b>BASKET PUNCH SLIM; UPBITER</b>						
	Upbiter straight	1.90	1.70	3.10	578-210-10	577-210-10
	Upbiter curved left 12°	1.90	1.70	3.10	578-220-13	577-220-13
	Upbiter curved right 12°	1.90	1.70	3.10	578-225-14	577-225-14
<b>BASKET PUNCH SCOOP STYLE</b>						
	Straight; small	2.10	1.60	3.20	578-310-00	577-310-00
	Straight; small; wide basket	2.10	1.60	4.00	578-335-00	577-335-00
	Straight; medium	2.30	1.80	4.10	578-350-00	577-350-00
	Straight; large; narrow basket	2.60	2.20	3.60	578-355-00	577-355-00
	Curved left 12°; small	2.10	1.60	3.20	578-365-03	577-365-03
	Curved right 12°; small	2.10	1.60	3.20	578-370-04	577-370-04
	Curved up 15°; small	2.10	1.60	3.20	578-320-02	577-320-02
	Curved down 15°; small	2.10	1.60	3.20	578-360-05	577-360-05
<b>BASKET PUNCH SCOOP STYLE; UPBITER</b>						
	Upbiter straight; small	2.10	1.70	3.20	578-315-10	577-315-10
	Upbiter straight; medium	2.30	2.20	4.10	578-340-10	577-340-10
	Upbiter curved left 12°; small	2.10	1.70	3.20	578-325-13	577-325-13
	Upbiter curved right 12°; small	2.10	1.70	3.20	578-330-14	577-330-14
	Upbiter curved up 15°; medium	2.30	2.20	4.10	578-345-02	577-345-02





- curved up  $15^\circ$
- curved down  $15^\circ$

curved right 12°

curved left 12°

straight



Description	Tip profile in mm	Tip width in mm	STELLA LOCK	
GRASPER ALLIGATOR; MAXIMUM				
Straight	2.60	3.90	578-720-00	
GRASPER ALLIGATOR; FINE				
Straight; fine	3.40	3.10	578-715-00	
GRASPER; 1 x 2 TEETH				
Straight; fine	4.00	3.90	578-710-00	
GRASPER; CUPPED				
Straight	3.30	2.80	578-730-00	

3



· WORLD WIDE DISTRIBUTION NETWORK ·

# KÖSLER

ARTHROSCOPY

**KOSLER** DSITRIBUTED BY OTHOTRUST MEDICAL

Tel: +1-844-9 878-7863  
1-844-9TRUSTME  
[orders@orthotrustmedical.com](mailto:orders@orthotrustmedical.com)

All rights are reserved. Technical alterations are possible. This folder may be used for no other purpose than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse, we retain the rights to recall our catalogs and pricelists and to take legal actions.