


neoligaments™



**Rapid Recovery**  
Orthopaedic Devices

Product Range

# Neoligaments

Neoligaments is a leading British medical device company specializing in the design, manufacture and supply of sterile implants for the sports medicine and orthopaedic markets.

The Neoligaments product portfolio comprises fixation devices and textile implants. The textile implants fall into two types: general devices, which are available in a range of standard sizes from which the surgeon can select the device best suited to repair the soft tissue injury; and speciality devices, which have been optimized to repair a particular tissue.


All our implants are single-use devices supplied within a double blister pack or double pouch pack and are sterilized by gamma irradiation.

[www.neoligaments.com](http://www.neoligaments.com)

## General Devices

### Poly-Tape

Poly-Tapes are intended to be used for tissue approximation, including use in reconstructing damaged or torn ligaments, tendons or other soft tissues.



### Ortho-Tape™

A high strength, narrow width, dense tape for tissue approximation, including use in reconstructing damaged or torn ligaments, tendons and other soft tissues. Also supplied for soft tissue (tendon and ligaments) fixation to bone with a Fastlok fixation device during orthopaedic reconstruction procedures.

- 102-1027** Ortho-Tape™ 3 x 600 mm (Multipack of 5)
- 102-1028** Ortho-Tape™ 3 x 400 mm (Multipack of 5)
- 102-1331** Ortho-Tape™ 3 x 400 mm, plus 8 mm Ti Fastlok™
- 102-1354** Ortho-Tape™ 3 x 400 mm (x2), plus 6 mm Ti Fastlok™



## Speciality Devices

### Rota-Lok™

For reconstruction of chronic, massive, full thickness rotator cuff tears.

102-1066



### AchilloCordPlus™

For repair of acute Achilles tendon ruptures.

102-1142



### JewelACL™

For reconstruction of the anterior cruciate ligament (ACL).

102-6003



### MPFL System™

For reconstruction of the medial patellofemoral ligament (MPFL). Implants include a 10 x 500 mm Poly-Tape and 6mm Ti Fastlok. Disposables include a probe and drill.

102-1060



## Fixation devices

For fixation of small tapes and sutures directly to bone.

### Fastlok™

For secure fatigue-resistant fixation of small tapes and multiple sutures directly to bone.

102-1380 6 mm Ti Fastlok

102-1381 8 mm Ti Fastlok



## Instruments

Fastlok instrument set includes:

### Impactor/Extractor

For insertion and removal of 6mm and 8mm Fastloks.

202-1137




### Sliding Hammer Assembly

To be used with the Impactor/Extractor for removal of 6mm and 8mm Fastloks.

202-1118





**Neoligaments™**  
A division of Xiros™  
Springfield House  
Whitehouse Lane  
Leeds LS19 7UE UK

Tel. +44(0)113 238 7202  
Fax. +44(0)113 238 7201  
[enquiries@neoligaments.com](mailto:enquiries@neoligaments.com)  
[www.neoligaments.com](http://www.neoligaments.com)

Registered in England No.1664824

Neoligaments, Xiros, Ortho-Tape,  
MPFL System, JewelACL, Rota-Lok,  
AchilloCord<sup>PLUS</sup> and Fastlok are  
trademarks of Xiros.

All rights reserved © Neoligaments™ 2010.  
Worldwide patents and patents pending.